

# Nanex™ multiSert+™

Models: NY1-SP/NC1-SP

## PERFORM MICRO INCISION CATARACT SURGERY WITHOUT COMPROMISES

### Nanex™ IOL

**Proven hydrophobic acrylic IOL material:** over 10 million units implanted worldwide, over 15 years<sup>1</sup>

**Significant PCO reduction:** proprietary active oxygen processing treatment<sup>2</sup>

### multiSert+™ injector

**Wound protection:**<sup>3</sup>

smallest nozzle size available for an open-loop preloaded hydrophobic IOL<sup>4</sup>

**Outstanding delivery consistency and smoothness:**

predictable lens release afforded by HOYA's unique patented preloaded IOL delivery mechanism<sup>5</sup>

**4-in-1 preloaded IOL system providing choice:**

push or screw injection, and *insert shield* in advanced or default position according to the surgeon's preference and individual case presented



The preloaded Nanex™ multiSert+™ system eliminates the compromises traditionally associated with sub-2.2 mm cataract surgery

# Nanex™ multiSert+™

Perform micro incision cataract surgery without compromises




## Nanex™ multiSert+™


Model name	NY1-SP	NC1-SP
<b>Optic design</b>	Aspheric design with sharp optic edge	
<b>Optic &amp; haptic material</b>	Hydrophobic acrylic with UV and blue light filter	Hydrophobic acrylic with UV-Filter
<b>Haptic design</b>	Modified C-loop, 5° angulation	
<b>Diameter (optic/OAL)</b>	6.0 mm/ 13.0 mm	
<b>Power</b>	+6.00 to +30.00 D (in 0.50 D increments)	
<b>Nominal A-constant*</b>	119.2	
<b>Injector</b>	multiSert+™ preloaded	
<b>Front injector tip outer diameter</b>	1.62 mm	
<b>Incision size</b>	As low as 1.8 mm	

\* The A-constant mentioned above is presented as a guideline only for lens power calculations. It is recommended that the A-constant measurement be customized based on the surgeon's experience and measuring equipment.



## multiSert+™



## Preloaded Nanex™ IOL



Nanex™ multiSert+™ is a 4-in-1 preloaded IOL system providing multiple choices at your fingertips



Information contained is intended for health care professionals. For a full list of indications and contra indications please refer to the Instructions For Use. Some of the products and/or specific features as well as the procedures featured in this document may not be approved in your country and thus may not be available there. Design and specifications are subject to change without prior notice as a result of ongoing technical development. Please contact our regional representative regarding individual availability in your country. HOYA, Nanex and multiSert+ are trademarks of the HOYA Corporation or its affiliates. ©2020 HOYA Medical Singapore Pte. Ltd. All rights reserved. HOYA Medical Singapore Pte. Ltd. | 455A Jalan Ahmad Ibrahim | Singapore 639939  
 HOYA Surgical Optics GmbH | De-Saint-Exupéry-Straße 10 | 60549 Frankfurt / Main | Germany Hotline DE: Tel. 0800 664 2 664 | Fax 0800 774 2 774 [hoyasurgicaloptics.com](http://hoyasurgicaloptics.com)

Singularly Focused. Globally Powered.™



0123

2020-09-30\_HSOE\_NY1-SP\_NC1-SP\_DS\_EN

**HOYA**  
SURGICAL OPTICS